

TargetAMD Symposium

**Quality Controls set-up during the
Development Phases of a Biological-
based Medicinal Product:**

WAS Gene Therapy Protocol

Salima Hacein-Bey-Abina

GTMP:“...active substance that contains or consists of a recombinant nucleic acid used in...view to regulate, repair, replace, add or delete a genetic sequence”.

Biological medicinal products:“... the active substance of which is a biological substance.

Biological substance:“...produced by or extracted from a biological source ..”.

Dir. 2001/83/EC , amd. by Dir. 2009/120/EC

Biologics: can be defined...by reference to their method of manufacture. *Guidelines for Good Manufacturing Practices for Medicinal Products for Human & Veterinary Use. Eudralex Volume 4*


ICH documents for Biologics

- **Q5 A:** Viral Safety
- **Q5 B:** Genetic Stability
- **Q5 C:** Product Stability
- **Q5 D:** Cell Substrates
- **Q5 E:** Comparability
- **Q6 B:** Specification
- **M4 / M2:** CTD / e-CTD
- **Q7:** GMP for APIs
- **Q8:** Pharmaceutical development
- **Q9:** Quality Risk Management
- **Q10:** Pharmaceutical quality system
- **Q11:** Development & Manufacture of Drug Substances


Wiskott-Aldrich Syndrome

- **Rare X-linked PID** (1-10/million births total)
- **Due to lack of WASP in blood cells**
 - cytoskeletal regulator via Arp2/ 3 complex activation
- **Complex symptoms**
 - Recurrent infections
 - Bleeding
 - Eczema
 - Auto-immune disorders
 - Increased incidence of malignancies
- **Treated by allogeneic BMT**


Genotypically id: 83%; related HLA-mismatched:50%



WASP Deficiency in the Immune System




- The WAS vector is an advanced-generation replication-defective self-inactivating lentiviral vector derived from rHIV1 and pseudotyped with VSVg.
- It encodes for the human WAS cDNA driven by its endogenous 1.6 kb promoter.



Clinical Trial for WAS: Phase I/II Pilot study

- Same vector and miscellaneous reagents used in three centers
- Severe WAS patients (score 3-5) with no HLA-geno-id BM or 10/10 or 9/10 Ag HLA-matched unrelated donor or cord blood unavailable
- Primary objectives:
 - safety and sustained engraftment of WASP⁺ transduced cells,
 - reconstitution of immunity
 - correction of thrombocytopenia
- Secondary objectives: to improve patient health (immunity and bleeding).





Applicable laws

Dir. 2001/20/EC, amend. Regulation No. 596/2009

Dir. 2005/28 /EC, Note for guidance on GCP,
CPMP/ICH/135/95, CPMP/ICH/291/95 - ICH E6&E8

Principles

Pre-requisite:

Benefits >risks, Patients understanding & consent,
Safeguards for physical & mental integrity, Insurance

Sponsor responsibility:

Trial conduct, investigator & institution, **QC (must be already set for both efficacy and safety evaluation before protocol starting)**, data collection stds, written protocol, investigator's brochure

GCP & Study design: ICH E6 & E8 Guidelines



Product characterization/Product Testing

Multiple preclinical steps must be defined and validated, the product should be characterized

- Define critical product attributes
- Define and control all components (cell types) present in the product
- Establish appropriate product specifications
 - Ensure the safety and consistency of product lots
 - Base acceptance criteria on experience
 - Agree with regulatory standards



Lentiviral -Transduced CD34+cells


Product characterization/Product Testing

Quality issues I

- **Non-clinical studies must allow to bring the first evidences of efficacy and safety through development of the most appropriate analytical tools adapted to the product under evaluation**
- **Starting material tested/documentated for viral, TSE safety**
- **Excipients, reagents and structural components qualified**
- **Product consistent and well characterised**
- **Active moiety identified and quantified**
 - **Manufacturing process validated and Potency assay available**
 - **Comparability issues**

Preclinical Safety data related to the selected WASP-LV vector

- Non-replicative, RCL-free
- Absence of vector toxicity *in vitro* on target cells
- Good safety profile from *in vivo* safety studies in mice
- No evidence of oncogenicity in IVIM assay
- Expected pattern of genomic integration



Modlich et al. Mol Ther 2009



Lentiviral -Transduced CD34+cells



Product characterization/Product Testing

Quality issues II

Control & validation of the manufacturing process

- **Critical steps & limits**
- Control via **relevant markers** (mRNA, protein expression..)

Purity


- Maximise active moiety
- Eliminate and control undesired components
- **Minimum requirement: consistency**

Tumourigenicity

- Limit amount of ...(dedifferentiated cells, **VCN**)
- Demonstrate genotypic/phenotypic stability during process

Gene Therapy Protocols Set-up

Quality issues



Viral batch tests: research & clinical grade

Viral titer, quality, potency

Pre-Runs:

***Production Process design, PTA**

Media, Cytokines, bags

***QC testing determination and validation**

***Specifications drawn up**

Target cells for GT:

Cell lines, CD34 CTRL, Patient CD34

Final corrected populations:

T-lymphocytes, DC, NK, B cells, granulocytes, Platelets

Hematopoietic potential of GM cells:

Clonogenic progenitors

LTC-IC

Efficacy of GT (Potency tests):

cell viability

cell number/dose

transduction rate

VCN (Q-PCR)

WASP Intra-cytoplasmic staining

Western Blot


In vitro prediction of efficacy

WASP restoration after Gene Transfer

Potency assays


- ✓ Quantitative or qualitative biological assay
- ✓ Predict safety/efficacy of the biologics under evaluation

Western blot analysis




Functional recovery (DC podosomes)
Restored Dendritic cells cytoskeletal dynamics


Intracellular immunostaining



WASP - differ. DC




Transduced diff. DC




Characteristics of Genetically Modified CD34+ cells

Patient FR01




CD34⁺ × 10⁶/kg	11
VCN in CD34	1.3

Patient FR02




CD34⁺ × 10⁶/kg	10.8
VCN in CD34	1.2

Mock TD2 (21%)



Mock TD1 (17%) TD2 (39%)



QC used to release to product

All the patients show a clear clinical benefit

Kinetics of mature WASP+ blood cells generation is tightly dependent on:

- the VCN**
- the number of CD34 infused**
- their own selective advantage**

Efforts are ongoing to improve and standardise the transduction procedure

Importance of reliable QC methods for product characterization since the clinical benefit is directly dependant on the « quality » of the genetically modified progenitor cells



Acknowledgements



M. Cavazzana

L. Caccavelli

J. Blondeau

C. Roudaut

B. Ternaux

C. Picard

A. Fischer

S. Blanche

N. Mahlaoui

A. Galy

S. Charrier

G. Honnet

F. Mavilio

A. Lim (I. Pasteur)

A. Thrasher, H. Gaspar (**London**)

A. Aiuti (**Milan**)

R. Bushman (**Philadelphia**)